

Hawkins Watts Australia Pty Ltd

Melbourne Office:
Suite 6, 2 Compark Circuit, Mulgrave Vic 3170
Phone 03 9561 3710

Sydney Office:
Suite 12, Level 1, 29 Kiora Road, Miranda NSW 2228
Phone 03 8543 3804

www.hawkinswatts.com.au



Carbonated Orange Drink using Instant CJ 204 Pectin

Recipe:

Instant CJ 204 Pectin (= 0.15 %)	1.5 g
Orange pulp (2x1)	10 g
Orange concentrate (6x1)	8.3 g
Orange flavour	0.5 g
Sodium benzoate	1 g
Artificial sweetener (e.g. Aspartame)	0.4 g
Citric acid, anhydrate	3.8 g
Water	28 g
Carbonated water	845 g

Input:	approx. 1000 g
Output:	approx. 1000 g
SS:	2 %
pH-Value:	3.4

Method:

- A. Mix all ingredients except carbonated water.
- B. Add A into water.
- C. Pasteurize.
- D. Fill
- E. Carbonate, close.